

## J-GEL



- ✓ I-gel is one of the new supraglottic airway tools.
- ✓ It is a device with a thermoplastic elastomeric structure designed not to press the laryngeal and pharyngeal anatomic structures.
- ✓ The cuff-like thickened structure has been developed from the mirror image of the supraglottic anatomy to adapt to that region.
- ✓ It has an additional lumen that allows aspiration of the stomach contents and a bite guard at the level of the guard.
- ✓ In this way, it is possible to evacuate the gas accumulation in the stomach during ventilation.
- ✓ Placement does not require much experience
- ✓ Because of a good fit in the larynx, oropharyngeal leak pressure is also high.
- ✓ There are seven different numbers that allow the use of children.
- ✓ Because it is not cuff, it is easy to place with its hard and flat structure, its placement time is short.
- ✓ There are also dimensions that allow the intubation tube to pass through.
- ✓ Tissue damage is less.